


LED Cold Light Source

LED3

**MADE
IN
GERMANY**

- Ideal replacement for 150W halogen cold light sources
- Passive, fanless cooling system
- Adaptable to various light guides
- Stackable and horizontal positioning
- Flicker-free lighting thanks to optimized electronic
- Very high luminous flux
- Brightness adjustable from 0 to 100%
- Solid aluminium housing in modern design



bestscientific.co.uk

For more information, contact Best Scientific
t: +44 (0)1793 812477 | e: enq@bestscientificweb.com

Technical Data

LED life span	50.000 h (according to manufacturer)
Light guide access	15 mm (Standard)
Power consumption	approx. 36 W
Input power supply	100 - 240 V AC, 50/60 Hz
Output power supply	15 V DC, 0 - 4 A
Dimensions (l × w × h)	259 mm × 174 mm × 81 mm
Material	anodized aluminium
Color	titanium and grey
Weight	approx. 3,7 kg

Versions

100-010055	LED3-P	passive cooling, 30 W, 5.700 K, CRI 70
100-010056	LED3-PC	passive cooling, 30 W, 5.700 K, CRI 90
100-010057	LED3-N	passive cooling, 30 W, 4.000 K, CRI 80
100-010058	LED3-W	passive cooling, 30 W, 3.000 K, CRI 80

Accessories

100-005989	Adaptor for Schott Light Guide	9 + 12 mm
100-006343	Adaptor for Schott Light Guide	10 + 12 mm

