



WI-FI CAMERAS

A unique solution to turn your old microscope into a Wi-Fi station, enabled to share your images with up to 6 devices at the same time. Simply logon with your tablet or smartphone and get ready to capture, edit and store your microscope images with the free App MotiConnect.

The Moticam X is equipped with its own on-board software, that can be easily controlled with a Web Browser which allows you to capture as well as measure (in case you are not allowed to download Apps). Alternatively, it also works with our standard Motic Images Plus Software for computers.

www.moticamseries.com | www.moticeurope.com

STREAMING WI-FI CAMERA



















Moticam X TECHNICAL SPECIFICATIONS

Sensor Type CMOS Sensor Size 1/3.2" 1.3 MP Resolution

Imaging Area 4.48mm x 3.36mm

Live Display Mode (through USB) 1280 x 960 (Wi-Fi only)

> Pixel Size 2.8µm x 2.8µm

Scan Mode Progressive

Wi-Fi **Data Transfer**

Exposure Time 0.05 msec to 0.07 sec

Operating Temperature From -10 to +60 Degrees Celsius non condensing

Max. Signal to Noise Ratio 42dB

> **Dynamic Range** 68dB

> > 1.0V/Lux-sec (550nm) Sensitivity

Support Device Universal Web browser (HTML5)

Supported OS Microsoft Windows XP/Vista/7/8/MAC OSX/Android 4,0 or higher/

iOS7 or higher

Minimum Computer Requirements 2GHz dualcore - RAM memory 2GB - Video memory min. 512 MB

> C-Mount Lens Mount Focusable Lens 12mm

Power supply 5V (USB - supply)

Package Includes CS Ring Adaptor, Focusable Lens, 30mm and 38mm Eyepiece

Adapters, USB 2.0 Cable, Calibration Slide, Macro Tube, Motic Images

Plus 3.0 for PC/OSX/Linux

*frames per second under optimal illumination conditions

MOTIC IMAGES PLUS 3.0 IMAGE ANALISYS SOFTWARE

Motic Images Plus 3.0

By connecting your Moticam to your computer you will be able to work with our well-known software, that comes included in the package. View, capture, edit, measure, make reports... all its standard features and the new ones packed in a new user-friendly interface.





MotiConnect App

View, edit, capture and store your images in the blink of an eye, view live images on multiple wireless devices and share your work with anyone in just a few clicks. The MotiConnect App makes data sharing possible in a collaborative environment, encouraging teamwork and collective learning experiences.

